PNEUMERICTM

CAPNOSPOTTM CE

INSTRUCTIONS FOR USE

Device Description

The Pneumatic Capnospot Pneumothorax Decompression Indicator is a capnography device that provides objective, real-time feedback during decompressive thoracostomy. Specifically, the device detects the presence of carbon dioxide in the gas from the decompressed pneumothorax, permitting more accurate placement of the needle or thoracostomy device compared to the current standard of care based upon auditory assessments.

The Pneumatic Capnospot is a therapeutic device intended to guide the placement of pneumothorax decompression needles and thoracostomy devices.

The product consists of a cylindrical polycarbonate tube containing a carbon dioxide detecting paper. The device has a male luer connector on the distal end of the device to connect to decompression needles and thoracostomy devices, a diaphragm to prevent backflow of gases, a color changing indicator and a female luer on the proximal end of the device that can be used to connect accessories.

The product is shipped non-sterile and labeled for single use only. The device is supplied in a single pack configuration labeled with a QR code linking to an electronic Instructions for Use. Each device is sealed in a foil pouch with a desiccant pack.


Figure 1. Pneumatic Capnospot Pneumothorax Decompression Indicator

Indication, Intended Use, and Expected Clinical Benefit

Used for more accurate placement of pneumothorax decompression devices than the current standard of care auditory assessments. The Pneumothorax Decompression Indicator is a colorimetric capnography device indicated to detect the carbon dioxide in the gas passing through the device used to treat pneumothorax with suspected tension physiology.

Contraindications

- Uncooperative patient.
- When pneumothorax decompression is otherwise contraindicated.

Warnings and Precautions:

- Read all instructions carefully. Failure to follow instructions may lead to serious medical consequences.
- Only persons having adequate training in pneumothorax decompression techniques should perform procedure.
- Capnospot has not been tested in patients with open pneumothorax, where treatment should be based on accepted guidelines and standard practices.
- Capnospot should be used in combination with clinical assessments and standard practice guidelines.
- Procedure details will vary according to practitioner's clinical judgement.
- Single Use Only.
- Do not use if package is opened or damaged.
- In order to achieve its intended use, Capnospot must be connected to a needle or thoracostomy device (angiocatheter, needle angiocatheter, pigtail thoracostomy catheter, or tube thoracostomy catheter with adapters) with a female luer fitting.
- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.
- Store between 10°F (-12°C) and 140°F (60°C). Storage for extended periods at temperatures below 10°F (-12°C) or above 140°F (60°C) may reduce shelf life.

Method for Use:

Prepare the Capnospot Pneumothorax Decompression Indicator for use:

1. Inspect the pouch and verify that it is unopened and undamaged. DO NOT use the device if pouch is found opened or damaged.
2. Gently open the pouch by tearing the notch below the seal and inspect the components for damage. DO NOT use the device if damage is evident.
3. The device consists of cylindrical polycarbonate tube with a male luer on the distal end, female luer on the proximal end, and blue indicator paper inside the tube.
4. A minimal amount of green-yellow color in the indicator paper is a normal condition. Do not use the device if majority of the indicator paper is yellow in color.
5. Connect the male luer on the Capnospot to the female luer portion of a needle or thoracostomy device as shown in figure 2 below.

NOTE: CAPNOSPOT MUST BE CONNECTED TO THORACOSTOMY DEVICE BEFORE PROVIDER ATTEMPTS TO DECOMPRESS PNEUMOTHORAX WITH SUSPECTED TENSION PHYSIOLOGY. ALTERNATIVELY, IF A NEEDLE ANGIOCATHETER IS BEING USED ONCE THE NEEDLE IS REMOVED, CAPNOSPOT SHOULD BE PROMPTLY ATTACHED TO THE ANGIOCATHETER.


Figure 2. Connection of Decompression Needle/Thoracostomy Device to Capnospot

6. Place the needle or thoracostomy device as clinically indicated for tension pneumothorax.
7. Observe for color change within 10 seconds. The indicator paper should change from blue to yellow if greater than ambient CO₂ levels are detected.
8. If color changes, continue with current clinical care.
9. If no color change occurs, keep the Capnospot connected to the needle or thoracostomy device, evaluate the clinical condition, and re-assess if needle decompression is indicated. Repeat steps 6-8 if required.

Important Risk Information

Potential adverse events associated with the use of Capnospot include:


- Delay of therapy

Disposal

- The carton is recyclable. Dispose of all packaging materials as appropriate.
- Use solid biohazard waste procedures to discard device.

Notice Regarding the Reporting of Adverse Events

Any serious incident that has occurred in relation to Capnospot should be reported to the manufacturer and the competent authority of the EU Member State or other relevant regulatory body in which the user and/or patient is established.

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Pneumeric, Inc.
OMC NI Medical Ltd
OMC Medical Limited
OMC Medical GmbH